


Today's Precision... Tomorrow's Vision

Latest Innovations from ASICO

AKAHOSHI MICRO-COAXIAL CAPSULORRHESIS FORCEPS WITH ROUND HANDLE AE-4344R




Angled 12.0mm jaws with special alignment mechanism allow for easy movement within the anterior chamber through incisions as small as 1.8mm without stretching the size. Triangular tips can be used as a cystotome. Round handle great for maneuverability. Laser lines 1.0mm apart.




TAN DALK SCISSORS – AE-5666/AE-5667

Safety platform with highly polished edges prevents perforation of Descemet's membrane. Available in right and left configuration.




ADVANCUT – CLEAR CORNEAL BLADE

The double bevel blade makes an ideal square incision. No need for hydration. The sharpest blade on the market. Available in multiple sizes with and without safety handle.



AKAHOSHI COMBO II PRECHOPPER – AE-4190

Reduces phaco time up to 50%. Can work with as small as 1.8mm incision. Sharp edge initiates incision and rounded blunt edge helps in safe division of nucleus at posterior capsule. Used for nuclei grades 1-2.



APAO BOOTH: 518 • ASCRS BOOTH: 2833-2825

The more you know...™

ASICO®

Today's Precision... Tomorrow's Vision®

ASICO LLC

26 Plaza Drive • Westmont, IL 60559 • USA

Phone: 630-986-8032 • Fax: 630-986-0065

Customer Service: 1-800-628-2879

www.asico.com • E-mail: info@asico.com

© 2011 ASICO LLC

AD86-11

CE