


Fixation Plates




Josephberg Trocar Fixation Plates

For Alcon's 23G and 25G trocar systems (Total Plus® Pak with Entry). Features two places for fixation. The U-shaped side-opening wraps around the collar of the trocar cannula to secure the cannula in place as the trocar is being removed. Available in two radii, 3.5mm (AE-2925C) and 4.0mm (AE-2926C). Sizes are laser-engraved on the top of the plate for confirmation under the microscope.




Total Plus® Pak is a trademark of Alcon Labs

Designed with Robert Josephberg, MD
Chief of Retina & Vitreous, New York Medical College
Valhalla, New York



Specialty Forceps



Tanaka Maculorrhexis Forceps: AE-4975



23G forceps with a specially designed tip to initiate the rhexis. Jaws of the forceps are treated for secure grasping.



Kawai ILM Forceps: AE-4980



23G forceps with a non-glare finish and specially designed platform for better visibility on the thin membrane in the macular region. Jaws are coated for secure grasping.



Kawai Sub Retinal Forceps: AE-4990



23G tapers to 25G at the tip. 4mm, serrated jaws help effectively grasp the sub retinal neo-vascular membrane.

Muscle Forceps




Jameson Muscle Forceps (right): AE-4410
Left configuration: AE-4411

Needle Holder




Barraquer Needle Holder (curved, with lock): AE-6160
Curved, without lock: AE-6161


Single-Use Retinal Cannulas




AS-7901
Packo Back Flush Needle
20G, blunt
Reusable Version: AE-7901



AS-7905
Packo Subretinal Fluid Drainage Needle
20G tapers to 25G, side opening
Reusable Version: AE-7905



AS-7903
Packo Back Flush Needle
20G tapers to 25G
Reusable Version: AE-7903




AS-7907
Subretinal Fluid Drainage Needle
20G, curved with silicone tip
Straight Version: AS-7906




Infusion Cannulas


AE-7960 2.5mm
AE-7961 4.0mm
AE-7962 6.0mm



Scleral Plugs



Comes in a package of 4



AE-7950
18G (black)
AE-7951
19G (red)

AE-7952
20G (gold)
AE-7953
Scleral Plug Set
(18, 19, and 20G)

Scleral Plug Forceps



Scleral Plug Forceps: AE-4900





Watzke Sleeve Spreader: AE-4910

Scleral Depressor



Packo Scleral Depressor: AE-1701



Flat handle on top represents marker position

Back Flush Handles



Packo Back Flush Handle: AE-7930
Manual back flush handle with control hole

Packo Back Flush Handle: AE-7931
Used with vitrectomy machine

The more you know...

ASICO

Today's Precision... Tomorrow's Vision

ASICO LLC
26 Plaza Drive • Westmont, IL 60559 • USA
Phone: 630-986-8032 • Fax: 630-986-0065
Customer Service: 1-800-628-2879
www.asico.com • E-mail: info@asico.com

© 2008 ASICO LLC

