


# Akahoshi Innovations

## Akahoshi Specula

- Small wire does not hinder the visibility of the surgical field or maneuverability of the instruments
- Adjustable mechanism easily secures the lid in place regardless of lid size
- 3-prong version is ideal to prevent cilia from escaping from the drape edge
- Short version has shorter prongs for deep set eyes and narrow lids



AE-1025  
AE-1026 (short)


## Akahoshi Sideport II Diamond Knife: AE-8104

- Trapezoidal blade
- 0.6mm at the shoulder, 1.0mm at the base



## Akahoshi Visco-Dissection Cannula: AE-7506

- Injects OVD between prechopped nuclei to separate from capsular bag
- Controls nuclear fragments during phaco
- Injects OVD after cataract has been removed
- Polishes the posterior capsule
- Provides counterforce during IOL implantation
- Injects OVD after IOL implantation to manipulate the lens in place



**ASICO**  
Today's Precision... Tomorrow's Vision™  
© 2009 ASICO LLC