

23G Vertical Scissors: AE-5609



Professor Shunji Kusaka
Department of Ophthalmology
Osaka University
Osaka, Japan

"These 23G vertical scissors work well with the 23G trocar system to dissect and peel fibrovascular membranes securely and safely. The curved blade and fine tip allow easy and atraumatic insertion into the space between the retina and fibrovascular membrane. The visibility and sharpness is excellent. Now with these scissors, I use the 23G trocar system in all of my diabetic retinopathy cases, even with extensive fibrovascular membranes."

Curved blade for excellent visibility!



- 45-degree angled scissors
- 2.3mm distal blade with thin profile allows precise control with minimal risk to choroids or underlying membrane

Featuring ASICO's micro V-R instrument design

Rounded handle for optimal tactile feedback

All-titanium body



Blue = forceps, spatula
Gold = scissors

Green: 23G
Pink: 25G
Silver: 27G

The more you know...™

ASICO®

Today's Precision... Tomorrow's Vision™

ASICO LLC
26 Plaza Drive • Westmont, IL 60559 • USA
Phone: 630-986-8032 • Fax: 630-986-0065
Customer Service: 1-800-628-2879
www.asico.com • E-mail: info@asico.com



© 2009 ASICO LLC

BR293-09