



# Espiritu-Caparas

## Chopper: AE-2574



**Cesar Ramon G. Espiritu, MD**  
**Manila, Philippines**

*"It made perfect sense for us to have both of our favorite 'second instruments' in one piece. We use the spatula to fixate the eye through the side port incision, manipulate and redirect lens material, and, by positioning it underneath the phaco tip, protect posterior capsule. The ideal curve, length and multiple cutting edges of the Verges-style chopper on the other end makes it a truly versatile one."*



- Verges-style horizontal chopper on one end: cuts through different nucleus densities on both sides as well as on its handle's axis and safely manipulates cataract fragments and intraocular structures
- 0.8mm spatula on other end: provides resistance against 'cheese-wiring' through softer lens material so one can move these pieces better and protect the posterior capsule by placing it underneath the phaco tip when emulsification is down to the last few pieces

*The more you know...*  
**ASICO**<sup>®</sup>  
*Today's Precision... Tomorrow's Vision*

ASICO LLC  
26 Plaza Drive • Westmont, IL 60559 • USA  
Phone: 630-986-8032 • Fax: 630-986-0065  
Customer Service: 1-800-628-2879  
www.asico.com • E-mail: info@asico.com

© 2009 ASICO LLC

