



2.2mm Capsulorhexis

Made Easy

No learning curve required! Our Kraff-Utratas have been modified for Micro Coaxial Cataract Surgery. With a much thinner profile, jaws open much wider without stretching the micro incision.

Thin Kraff-Utrata Capsulorhexis Forceps: AE-4503RH



- Sharp tips can be used as cystotome
- Jaws are recessed to prevent iris pickup
- 1.0mm laser lines aid in measuring capsulorhexis




“These forceps are very easy to use through a 2.2mm incision.”

Sheridan Lam, MD, Elmhurst, Illinois

Thin Vaulted Kraff-Utrata Capsulorhexis Forceps: AE-4504RH



- Sharp tips effectively grasp the anterior capsule
- Jaws are vaulted to conform to the curvature of the lens
- 1.0mm laser lines aid in measuring capsulorhexis



The more you know...

ASICO

Today's Precision...Tomorrow's Vision

ASICO LLC
26 Plaza Drive • Westmont, IL 60559 • USA
Phone: 630-986-8032 • Fax: 630-986-0065
Customer Service: 1-800-628-2879
www.asico.com • E-mail: info@asico.com

© 2009 ASICO LLC



BR299-09