


OR Essentials

Castroviejo Caliper, straight:
AE-1500

Overall length: 86mm
Measures 1.0-20.0mm
in 1.0 increments



Castroviejo Lacrimal Dilator:
AE-1640

Overall length: 140mm



Infant Lacrimal Dilator:
AE-1650

Overall length: 70mm
Stainless steel



May Scleral Depressor:
AE-1702

Overall length: 93mm
Sandblasted tip



Fine-Thorton Fixation Ring:
AE-2198

Overall length: 115mm
13.0mm diameter with 8.0mm opening and atraumatic fixation teeth



Sinskey II IOL Manipulating Hook:
AE-2208

Overall length: 115mm; Tip diameter: 0.2mm
Blunt, angled tip



Kuglen Iris Hook & IOL Manipulator:
AE-2231

Overall length: 120mm; Tip diameter: 0.80mm
Push-pull model with angled tip



Maloney No Hole IOL Manipulator:
AE-2413

Overall length: 120mm; Tip diameter: 0.25mm
Semi-angled tip



Algerbrush Rust Ring Remover, 0.5mm:
AE-3250

Complete set with
0.5mm burr
Replacement Burr:
AE-3251




Harms-Colibri Corneal Forceps
AE-4036

Overall length: 70mm
Very delicate, 0.12mm tips with 1x2 teeth and
6mm tying platform




Bishop-Harm on Tissue Forceps, delicate:
AE-4101

Overall length: 86mm
Delicate, 0.8mm tips with 1x2 teeth



Bishop-Harm on Dressing Forceps:
AE-4102

Overall length: 85mm
Criss-cross serrated tips



Harm's Tying Forceps, straight:
AE-4350

Overall length: 102mm
Smooth jaws with 6.0mm tying platform and 0.5mm tips
Ideal for 7-0 and 9-0 sutures



Jaffe Tying Forceps, straight:
AE-4360

Overall length: 102mm
Smooth jaws with 6.0mm tying platform and 0.3mm tips
Ideal for 8-0 and 11-0 sutures



Jaffe Tying Forceps, curved:
AE-4361

Overall length: 102mm
Smooth jaws with 6.0mm tying platform and 0.3mm tips
Ideal for 8-0 and 11-0 sutures




McPherson Tying Forceps, straight:
AE-4370

Overall length: 82mm
Smooth jaws with 5.0mm tying platform and 0.3mm tips
Ideal for 8-0 and 11-0 sutures




McPherson Tying Forceps, angled:
AE-4371

Overall length: 82mm
Smooth jaws with 5.0mm tying platform and 0.3mm tips
Ideal for 8-0 and 11-0 sutures




Kelman-McPherson Forceps:
AE-4372

Overall length: 82mm
Smooth jaws with 7.5mm tying platform and 0.3mm tips
Ideal for implanting IOLs



McPherson Tying Forceps, angled, long handle:
AE-4374

Overall length: 102mm
Smooth jaws with 5.0mm tying platform and 0.3mm tips
Ideal for 8-0 and 11-0 sutures



ASICO Tying Forceps, straight:
AE-4377

Overall length: 88mm
5.0mm tying platform and precise finish on tips



Tennant Tying Forceps, straight:
AE-4380

Overall length: 102mm
Smooth jaws with 6.0mm tying platform and 0.25mm tips
Ideal for 9-0 and 11-0 sutures




Tennant Tying Forceps, curved:
AE-4381

Overall length: 102mm
Smooth jaws with 6.0mm tying platform and 0.25mm tips
Ideal for 9-0 and 11-0 sutures



Gaskin Fragment Forceps:
AE-4391

Overall length: 102mm
Angled, very delicate, smooth jaws



Berke Ptosis Forceps, 20mm:
AE-4440

Overall length: 107mm
Longitudinally-grooved jaws and sliding lock



Barraquer Glia Forceps:
AE-4520

Overall length: 103mm
Slightly angled, smooth jaws



Castroviejo Universal Corneal Scissors, small
AE-5052

Overall length: 108mm
Curved, 10.0mm blades
with blunt tips



Stern-Gills Capsulotomy Scissors:
AE-5490

Overall length: 108mm
Extra thin, 10.0mm blades



McPherson-Westcott Conjunctival Scissors:
AE-5505

Overall length: 108mm
Curved, 10.0mm blades
with blunt tips



McPherson-Westcott Stitch Scissors:
AE-5509

Overall length: 108mm
Very sharp, 10.0mm blades with pointed tips



Barraquer Needle Holder, curved with lock:
AE-6140

Overall length: 135mm
Delicate, 9.0mm jaws
Ideal for 8-0 to 11-0 sutures




Castroviejo Needle Holder:
AE-6220

Overall length: 140mm
Light serrated, 10.0mm jaws
Ideal for 5-0 to 6-0 sutures




Rycroft Air Injection Cannula:
AE-7000 (30G) / AE-7001 (27G)

Overall length: 32mm
Angle to tip: 5.0mm




Lacrim al Cannula, 23G, curved:
AE-7011

Overall length: 32mm
12mm malleable tip




Atkinson Retrobulbar Needle, 25G:
AE-7071

Overall length: 35mm
Beveled tip



Bishop-Harm on Irrigating Cannula, 19G:
AE-7200

Overall length: 34mm
Angle to tip: 7.0mm



ASICO Double Male Adapter, Teflon®
AE-7389

Overall length: 50mm



Silicone Tubing, 12in:
AE-7332 (large) / AE-7333 (small)

Large: 0.078x0.024
Small: 0.030x0.018

