



Inamura II Hydrodissection Cannula



Mikio Inamura MD
Yokohama, Japan



"With this new cannula I designed with ASICO, there is no difficulty in performing hydrodissection. I recommend using the cannula with a 5 or 10cc syringe. Softly insert the tip right beneath the edge of the CCC and inject BSS. Make sure that the tip is not too far from the edge of the CCC."



DVD
AVAILABLE



Available in 22G (AE-7642) and 25G (AE-7644)

Advantages:

- 45° angled flow for efficient hydrodissection
- Sandblasted tip increases control in the presence of visco especially in cases where the nucleus is rotated inside the capsular bag.

The more you know...
ASICO[®]
Today's Precision...Tomorrow's Vision™

ASICO LLC
26 Plaza Drive • Westmont, IL 60559 • USA
Phone: 630-986-8032 • Fax: 630-986-0065
Customer Service: 1-800-628-2879
www.asico.com • E-mail: info@asico.com
© 2009 ASICO LLC

