


Phaco Instruments

Specula


Barraquer Wire Speculum, large
AE-1002

Overall length: 40mm
Blade length: 14mm



Barraquer Heavy Wire Speculum
AE-1003

Overall length: 40mm
Blade length: 14mm
1.0mm heavy wire



Miyoshi Adjustable Small Orbit Speculum
AE-1014N

Overall length: 70mm
Stainless steel



Capsulorhexis Forceps

Inamura Super III Capsulorhexis Forceps
AE-4346



Patented reverse cross-action mechanism
minimizes incision stretch and reduces visco loss



Thin Vaulted Kraff-Utrata Capsulorhexis Forceps
AE-4504RH


Thin profile jaws maneuver easily within a small
incision and reduces visco loss



Blades


Kellan Disposable Blades, 2.5mm
AS-8010AH

2.5mm Angled Blade
6/box



Kellan Disposable Blades, 3.2mm
AS-8011AH

3.2mm Angled Blade
6/box



Hydrodissection Cannulas

Nagahara Hydrodissection Cannula II:
AE-7622A

27G

Flat tip allows targeted dispersion of BSS for
effective hydrodissection




Miyoshi Vertical Ultimate Dual-Flow
Hydrodissection Cannula: AE-7643



23G

120° angled flow allows complete
hydrodissection the very first time



Oki/Yoshitomi Hydrodissection Cannula
AE-7648

23G

Precisely penetrates between the capsular bag
and the epinucleus for complete hydrodissection



Prechoppers

Akahoshi Combo II Prechopper
AE-4190



designed for nuclei grades 1-2
sharp side for initial crack
blunt side to safely complete division



Akahoshi Universal II Prechopper
AE-4192



designed for hard nuclei with the
Counter Prechop Technique
sharp tip for easy penetration



Hill Prechopper
AE-4193



designed for nuclei grades 2-4
hook design for easier initial crack
blunt side to safely complete division




Choppers

Koch Spatula/Nagahara Chopper
AE-2520



spatula used to manipulate and lift
chopper horizontally chops toward
phaco handpiece



Koch-Minami Chopper
AE-2528

versatile instrument acts as a divider, dissector,
manipulator, and protector



Nagahara Karate Chopper
AE-2529

spear tip with inner cutting edge vertically
chops and pushes nucleus towards
phaco handpiece



I/A Handpieces, Cannulas, and Tips

Micro Irrigating Handpiece, curved
AE7-0208



23G
dual 0.4mm oval ports



Micro Aspirating Handpiece
AE7-0217/AE7-0218 (w/polisher)




23G
0.3mm port



Viscoelastic Aspirating Cannula
AE-7270

22G
Overall length: 32mm
Angle to tip: 10mm



Akahoshi Mini Straight I/A Tip
AE7-3050



22G
small 0.70mm diameter allows ample flow for
superior chamber maintenance



The more you know...™
ASICO
Today's Precision... Tomorrow's Vision

ASICO LLC
26 Plaza Drive • Westmont, IL 60559 • USA
Phone: 630-986-8032 • Fax: 630-986-0065
Customer Service: 1-800-628-2879
www.asico.com • E-mail: info@asico.com
© 2009 ASICO LLC



BR333-09