



Bimanual I/A Handpieces

NEW!

22 Gauge




Irrigation: AE7-0110

- Laser welded tip gives better durability during sterilization
- Ergonomic handle for comfort
- Spatulated tip allows easy access
- Front opening allows maximum flow for better AC maintenance



Aspiration: AE7-0120

- Laser welded tip gives better durability during sterilization
- Ergonomic handle for comfort
- Curved with capsule polisher
- 0.3mm port




23 Gauge





Irrigation: AE7-0207 (Straight) AE7-0208 (Curved)

- Laser welded tip gives better durability during sterilization
- Ergonomic handle for comfort
- Dual oval ports 0.3mm x 0.4mm



Aspiration: AE7-0217 (Curved) AE7-0218 (Curved w/Capsule Polisher)

- Laser welded tip gives better durability during sterilization
- Ergonomic handle for comfort
- 0.3mm port



*The more you know...*TM
ASICO[®]
 Today's Precision... Tomorrow's Vision[®]

ASICO LLC
 26 Plaza Drive • Westmont, IL 60559 • USA
 Phone: 630-986-8032 • Fax: 630-986-0065
 Customer Service: 1-800-628-2879
 www.asico.com • E-mail: info@asico.com

