


ASICO I/A Handpiece


AE7-0029

The only I/A handpiece in the market to disassemble for thorough cleaning

Handle comes completely apart for easier cleaning



Universal threading allows to use this handle with various threaded I/A tips



ASICO I/A Handpiece:

- Compatible to all popular phaco consoles
- Designed to be used with various I/A tips
- Fits multiple sleeves
- All titanium

The more you know...
ASICO[®]
Today's Precision... Tomorrow's Vision

ASICO LLC
26 Plaza Drive • Westmont, IL 60559 • USA
Phone: 630-986-8032 • Fax: 630-986-0065
Customer Service: 1-800-628-2879
www.asico.com • E-mail: info@asico.com
© 2011 ASICO LLC BR 436-11

CE