


Your Choice of Toric Markers

Prior to Cataract Surgery...




Nuijts/Lane Pre-op Toric Reference Marker: AE-2791TBL

Bubble-level for precise alignment
Used temporarily



Pre-op Toric Reference Marker: AE-2793D

Lightweight, titanium body
Used temporarily



Pre-op Toric Reference Marker, Single-Sided: AE-2793S

Handle position allows even marking with
less pressure



Prior to Toric IOL Implantation...



Nuijts/Solomon Toric Axis Marker: AE-2740N

Circular opening for maximum visibility
Used with the Beveled Degree Gauge



Maloney Astigmatism Axis Marker: AE-2741

Laser lines indicate marking pattern
Used with the Beveled Degree Gauge



ASICO Beveled Degree Gauge: AE-1590

Beveled design for maximum visibility
under the microscope



Modi Intra-Op Toric Marker: AE-2848

Designed for small eyes
Very thin profile



Intra-Op Toric Axis Marker II: AE-2794

Marker and degree gauge in one
Cutouts for precise alignment



The more you know...™
ASICO
Today's Precision... Tomorrow's Vision

ASICO LLC
26 Plaza Drive • Westmont, IL 60559 • USA
Phone: 630-986-8032 • Fax: 630-986-0065
Customer Service: 1-800-628-2879
www.asico.com • E-mail: info@asico.com
© 2008 ASICO LLC

