

Today's Precision...Tomorrow's Vision

Introducing ICMS


Intraocular Complication Management Set



Alan Crandall, MD
Salt Lake City
Utah, USA

"Every center should have a set on hand to deal with complicated cases. Standard instruments may be used 'reasonably well' but when I encounter cases with complications, I want to offer my patients the best."

The ICMS instrument set is designed by Dr. Alan Crandall of Salt Lake City, Utah, USA.



AE-6902 ICMS 21 Gauge Needle Holder Straight

Specially treated tungsten carbide jaws create a secure hold on the needle allowing for an easy Intraocular suturing procedure.



AE-6903 ICMS 21 Gauge Needle Holder Curved

Specially treated tungsten carbide jaws create a secure hold on the needle allowing for an easy Intraocular suturing procedure.



AE-4906 ICMS 23 Gauge Tying Forceps Straight

AE-4908 ICMS 23 Gauge Tying Forceps Curved

Features specially designed jaws for handling the sutures securely without damage.



AE-5613 ICMS 21 Gauge Scissors Curved

A multi-functional design that features curved blades to effectively cut capsular tags, irises, haptics of IOLs and sutures.



These instruments feature an ergonomic color-coded easy to identify titanium handle with an easy to use protective sterilizable cover. The collection includes 21 & 23 gauge state of the art designs.

The more you know...™

ASICO®

Today's Precision...Tomorrow's Vision™

ASICO LLC
26 Plaza Drive • Westmont, IL 60559 • USA
Phone: 630-986-8032 • Fax: 630-986-0065
Customer Service: 1-800-628-2879
www.asico.com • E-mail: info@asico.com

© 2008 ASICO LLC



BR182-08