



José Güell, MD
Barcelona, Spain




Dimensions

Overall length: 102mm
Tip-to-bend: 9.5mm
Width at tip: 1.0mm

AE-4210: Güell DMEK Forceps

"The Güell DMEK forceps are highly polished and offer a broad grasping surface to peel the endothelium membrane with one hand, diminishing the risk of tearing the membrane. Of all the forceps I have tried for endothelium dissection, these are by far the safest and most effective."

Güell DMEK Forceps in Action



Available on...



[youtube.com/ASICOLLC](https://www.youtube.com/ASICOLLC)

Scan to Watch Video



The more you know...
ASICO®
Today's Precision... Tomorrow's Vision™

ASICO LLC

26 Plaza Drive • Westmont, IL 60559 • USA

Phone: 630-986-8032 • Fax: 630-986-0065

Customer Service: 1-800-628-2879

www.asico.com • E-mail: info@asico.com

©2012 ASICO LLC

