

Modi Femtosecond Speculum

Dr. Modi's NEW INVENTION




Dr. Satish Modi
Fishkill, NY

Silastic tubing on the tip



AE-1037
Modi Femtosecond Speculum

Silastic tubing on the tip



AE-1057
Modi Femtosecond Speculum, V-style

PRODUCT ADVANTAGES

- Ergonomically perfect for a dock on each and every patient
- What makes it unique are its silastic coated open blades
 - they are adjustable
 - they allow lax conjunctiva to be kept gently away, to allow a dock, even in patients with floppy lids
 - they allow backwards pressure on the globe, to enable a dock in patients with easily ballotable eyes
 - they ensure that the conjunctiva doesn't get dragged into the femtosecond cone